Test report

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Sample Information

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Sample Name 样品名称	Norepinephrine Bitartrate Injection 重酒石酸去甲肾上腺素注射液	Sample properties 样品性状	Liquid 液体
Date of receipt 收样日期	2025/10/30	During testing 测试期间	2025/10/30 ~ 2025/11/05

Test the ingredients and structural formula

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Experimental Requirements

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Under the premise that the detection method parameters and chromatographic column specifications (4.6 × 250 mm, 5 μm) remain unchanged, screen for other suitable C18 chromatographic columns. The specific requirements are as follows:

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The chromatographic peak at around 2.5 minutes (i.e., the solvent peak) shall be a sharp peak and have good resolution from the chromatographic peak at around 3 minutes.

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The main peak shall elute at approximately 6 minutes (±1 minute).

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The selected chromatographic columns shall be easy to use and store, and columns requiring refrigerated storage shall be excluded.

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Tel:400-810-6969

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Web:www.welchmat.com

Chromatographic column: 色谱柱	十八烷基硅烷键合硅胶 (4.6*250mm, 5 μ m)		
Mobile phase A: 流动相 A:	<p>Take approximately 15.601 g of sodium dihydrogen phosphate (ranging from 14.0409 g to 17.1611 g) and about 0.5 g of sodium heptanesulfonate (ranging from 0.45 g to 0.55 g). Dissolve them in water and dilute the solution to 1000 milliliters. Adjust the pH value to 3.0\pm0.2 with phosphoric acid. After filtering through a 0.45 μm aqueous filter membrane and sonicating for 10 minutes, the desired solution is obtained.</p> <p>取磷酸二氢钠约 15.601g (14.0409g~17.1611g), 庚烷磺酸钠约 0.5g (0.45g~0.55g), 加水溶解至 1000ml, 用磷酸调节 pH 值至 3.0\pm0.2, 经 0.45μm 水膜滤过后, 超声 10 分钟后, 即得。</p>		
Mobile phase B: 流动相 B:	<p>Take 990 milliliters of methanol, add 10 milliliters of isopropanol, mix well. After filtering through a 0.45 μm organic filter membrane and sonicating for 10 minutes, the desired solution is obtained.</p> <p>取甲醇 990ml, 加异丙醇 10ml, 混匀, 经 0.45μm 有机滤膜过滤后, 超声 10 分钟后, 即得。</p>		
Flow rate: 流速	1.0 mL/min		
Injection volume: 进样量	20 μ L		
Column temperature: 柱温	25 $^{\circ}$ C		
Detector: 检测器	UV		
Detection wavelength: 检测器波长	280nm		
Gradient program 梯度程序	Gradient Elution 梯度洗脱		
	T (min)	流动相 A	流动相 B
	0	100	0
	4	100	0
	10	99	1
	20	70	30
	25	50	50
	25.1	100	0
35	100	0	
Reagent information 试剂信息			
Reagent Name 试剂名称	Grade 级别	Brand 品牌	
Water 水	Ultrapure water 超纯水	Welch 月旭	
Methanol 甲醇	HPLC 色谱级	Welch 月旭	

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1.2. Solution preparation (溶液配制)**1.2.1. Preparation of mobile phase (流动相配制):**

Mobile phase A: Take approximately 15.601 g of sodium dihydrogen phosphate (ranging from 14.0409 g to 17.1611 g) and about 0.5 g of sodium heptanesulfonate (ranging from 0.45 g to 0.55 g). Dissolve them in water and dilute the solution to 1000 milliliters. Adjust the pH value to 3.0±0.2 with phosphoric acid. After filtering through a 0.45 μm aqueous filter membrane and sonicating for 10 minutes, the desired solution is obtained.

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Mobile phase B: Take 990 milliliters of methanol, add 10 milliliters of isopropanol, mix well. After filtering through a $0.45 \mu\text{m}$ organic filter membrane and sonicating for 10 minutes, the desired solution is obtained.

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None

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
1.2.3. Test solution (测试溶液):

Take this solution and inject directly.

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2. Chromatogram and data (谱图和数据)

(1) Analysis spectrum of Test solution (测试溶液分析图谱):

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Tel: 400-810-6969

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Web: www.welchmat.com

Figure 1 Typical Chromatogram of Single Injection
 图 1 单个进样典型图


Figure 2 Comparison Chart of Test Results from Different Batches
 图 2 不同批号测试结果对比图


Figure 3 Comparison Chart of Repeatability Test Results (n=6)

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Tel:400-810-6969

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Web:www.welchmat.com


Figure 4 Collection Results (Taking a Single Collection Result as an Example)

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3. Conclusion (结论)

When the Yuexu chromatographic column Ultimate AQ - C18 (4.6 × 250 mm, 5 μm) is used under these chromatographic conditions, the main peak elutes with a retention time of approximately 6 minutes (±1 minute). The solvent peak presents a sharp shape, and its resolution from the impurity peak at around 3 minutes meets the requirements. Meanwhile, this chromatographic column is easy to store and requires no refrigeration. A repeatability test was conducted using this column, and the test results indicated good repeatability. In conclusion, this chromatographic column meets all the aforementioned requirements of the customer.

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