Test report

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| Sample Information   |                               |                           |                       |
|--|-------------------------------|---------------------------|-----------------------|
| 样品信息   |                               |                           |                       |
| Sample Name<br>样品名称  | Astragalus membranaceus<br>黄芪 | Sample properties<br>样品性状 | Liquid<br>液体          |
| Date of receipt<br>收样日期  | 2025/08/19                    | During testing<br>测试期间    | 2025/09/18-2025/09/21 |
| Test the ingredients and structural formula  |                               |                           |                       |
| 测试成分及结构式   |                               |                           |                       |
| Astragaloside IV<br>黄芪甲苷   |                               |                           |                       |
| Experimental Requirements  |                               |                           |                       |
| 实验要求   |                               |                           |                       |
| The resolution in system suitability solution A shall not be less than 2.0 (calculated as the ratio of the polymer peak height to the peak-to-trough height between liraglutide and the polymer)..<br>按客户方法筛柱子, 要求分离度 1.5 以上 |                               |                           |                       |
| Reference Method   |                               |                           |                       |
| 参考方法   |                               |                           |                       |
| General Chapters (Volume IV) 0512 of the Chinese Pharmacopoeia 2025 Edition<br>中国药典 2025 年版四部通则 0512   |                               |                           |                       |
| Reagent information  |                               |                           |                       |
| 试剂信息   |                               |                           |                       |
| Reagent Name<br>试剂名称   | Grade<br>级别                   | Brand<br>品牌               |                       |
| Water<br>水   | Ultrapure water<br>超纯水        | Welch<br>月旭               |                       |
| Glacial acetic acid<br>乙腈  | HPLC<br>色谱级                   | Welch<br>月旭               |                       |

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Both the test solution of Astragalus membranaceus (Huangqi) and the reference standard solution of Astragaloside IV are provided by the customer, and can be directly used for sample injection analysis.

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

### (6) Tset Solution Full Range Profile (供试品溶液全量程图谱)

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
### 3. Conclusion (结论)

Using Welch Ultimate Plus C18 (4.6×250mm,5μm), the detection requirements can be met under these chromatographic conditions.

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