## Test report

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Sample Information			
样品信息			
Sample Name 样品名称	Semaglutide 司美格鲁肽	Sample properties 样品性状	Solid 固体
Date of receipt 收样日期	2025/09/29	During testing 测试期间	2025/10/20 ~ 2025/10/22
Test the ingredients and structural formula			
测试成分及结构式			
API/主成分: BOC -His(Trt)- Aib- Glu(OtBu)-GlyOH Isomer/异构体: BOC-D-His(Trt) -Aib-D- Glu(OtBu)-GlyOH			
Experimental Requirements			
实验要求			
Develop a method for determining the content of polypeptides, requiring the separation of the main component from its isomers with a resolution greater than 2. 开发方法测定多肽含量, 要求主成分与异构体分离, 分离度大于 2			
Reference Method			
参考方法			
/			
Reagent information			
试剂信息			
Reagent Name 试剂名称	Grade 级别	Brand 品牌	
Water 水	Ultrapure water 超纯水	Welch 月旭	
Acetonitrile 乙腈	HPLC 色谱级	Welch 月旭	
Trifluoroacetic acid	AR	Aladdin	

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#### 1.2. Solution preparation (溶液配制)

##### 1.2.1. Preparation of mobile phase (流动相配制):

Mobile phase A :Take 1ml TFA in 1000 ml of pure water, mix well and filter to obtain.

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Mobile phase B: Take 1ml TFA in 1000 ml of acetonitrile, mix well and filter to obtain.

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1.2.2. Blank solution、Diluent (空白溶液、稀释剂):

Acetonitrile-Water (1:1, v/v) (等比例乙腈水)

1.2.3. API reference solution: Weigh 2 mg of the main component, transfer it to a 10 mL volumetric flask, dissolve it with a solvent, dilute to the marked volume, and mix well to obtain the solution.

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1.2.4. Isomer reference solution: Weigh 1 mg of the isomer, transfer it to a 10 mL volumetric flask, dissolve it with a solvent, dilute to the marked volume, and mix well to obtain the solution.

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1.2.5. Mixed reference solution: Take the reference solutions of the main component and the isomer, mix them in equal proportions, and blend well to obtain the solution.

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
1.2.6. System solution: Weigh 10 mg of the main component accurately, transfer it to a 10 mL volumetric flask, dissolve it with a solvent, add 0.1 mL of the isomer reference solution, dilute to the marked volume, and mix well. This gives a system solution with a main component concentration of 1 mg/mL and an isomer concentration of 1  $\mu$ g/mL.

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
## 2. Chromatogram and data (谱图和数据)

(1) Analysis spectrum of blank solution (空白溶液分析图谱):


(2) Analysis spectrum of Sample API reference solution (主成分对照溶液分析图谱):

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(3) Analysis spectrum of Isomer reference solution (异构体对照溶液分析图谱):

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
#### (4) Analysis spectrum of Mixed reference solution (混合对照溶液分析图谱):

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
#### (5) Analysis spectrum of System solution (系统溶液分析图谱)

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
### 3. Conclusion (结论)

Using Welch Ultimate LP-C18 (4.6×250mm, 5µm) analyzed the content of polypeptide under this chromatographic condition, the resolution between the main component and the isomer is 6.1, which was analyzed without interference from other substances, meets the customer's testing requirements.

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