Test report

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Sample Information			
样品信息			
Sample Name 样品名称	Test Sample 供试品	Sample properties 样品性状	Solid 固体
Date of receipt 收样日期	2025/10/20	During testing 测试期间	2025/10/23 ~ 2025/10/28
Test the ingredients and structural formula			
测试成分及结构式			
GHB-7-30		GHB-8-3	
		ZMC-10-52	
Experimental Requirements			
实验要求			
<p>Develop a method to analyze the purity of ZMC-10-52 in the test sample, as well as whether GHB-7-30 and GHB-8-3 remain in it.</p> <p>开发方法分析供试品中 ZMC-10-52 的纯度, 以及是否残留 GHB-7-30 及 GHB-8-3</p>			
Reference Method			
参考方法			
/			
Reagent information			
试剂信息			
Reagent Name 试剂名称	Grade 级别	Brand 品牌	
Water 水	Ultrapure water 超纯水	Welch 月旭	

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## 1.2. Solution preparation (溶液配制)

### 1.2.1. Preparation of mobile phase (流动相配制):

Mobile phase A: Measure 0.8 mL of phosphoric acid into 1000 mL of purified water, mix thoroughly, and then perform suction filtration to obtain the desired solution.

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Mobile phase B: Obtain the product by suction filtration with Acetonitrile.

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### 1.2.2. Blank solution、Diluent (空白溶液、稀释剂):

Acetonitrile-Water (1:1, v/v) (等比例乙腈水)

1.2.3. Test solution: Weigh 5 mg of the main component, transfer it to a 5 mL volumetric flask, dissolve it with a solvent, dilute to the marked volume, and mix well to obtain the solution.

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1.2.4. GHB-7-30 reference solution: Weigh 5 mg of the reference substance, transfer it to a 5 mL volumetric flask, dissolve it with a solvent, dilute to the marked volume, and mix well to obtain the solution.

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1.2.5. GHB-8-3 reference solution: Weigh 5 mg of the reference substance, transfer it to a 5 mL volumetric flask, dissolve it with a solvent, dilute to the marked volume, and mix well to obtain the solution.

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


(3) Analysis spectrum of GHB-8-3 reference solution (GHB-8-3 对照溶液分析图谱)

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Tel:400-810-6969

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
Web:www.welchmat.com

(4) Analysis spectrum of Test solution (供试品溶液分析图谱)

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Web: www.welchmat.com

### 3. Conclusion (结论)

Using Welch Welchrom C18 (4.6×250mm, 5μm) analyzed test substances under this chromatographic condition, the purity of its main component is 83.3%, which meets the customer's testing requirements.

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