Test report

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Sample Information			
样品信息			
Sample Name 样品名称	Norfloxacin Lactate 乳酸诺氟沙星	Sample properties 样品性状	Solid 固体
Date of receipt 收样日期	2025/09/11	During testing 测试期间	2025/09/12 ~ 2025/09/18
Test the ingredients and structural formula			
测试成分及结构式			
 Norfloxacin Lactate 乳酸诺氟沙星			
Experimental Requirements			
实验要求			
<p>The peak shape of norfloxacin needs to be improved. The blank solvent shows no interference. In the chromatogram of the reference solution, the theoretical plate number of the main peak, calculated based on norfloxacin, should be not less than 1500.</p> <p>诺氟沙星峰峰型需改善, 空白溶剂无干扰, 对照品溶液图中, 主峰理论塔板数按诺氟沙星计算应不低于 1500。</p>			
Reference Method			
参考方法			
/			
Reagent information			
试剂信息			
Reagent Name	Grade	Brand	

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Web: www.welchmat.com

1.2.1. Preparation of mobile phase (流动相配制):

0.025 M Phosphoric Acid Solution: Weigh 2.45 g of phosphoric acid, dissolve it in pure water to a total volume of 1000 ml, and after dissolving thoroughly, adjust the pH to 2.2 with triethylamine; the solution is obtained.

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Measure 870 ml of 0.025 M phosphoric acid solution and 130 ml of ultrapure water, mix them thoroughly, then perform suction filtration to obtain the solution.

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1.2.2. Blank solution、Diluent (空白溶液、稀释剂):

Pure water.

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0.1 M Hydrochloric Acid Solution: Measure 9 ml of hydrochloric acid, place it in a 1000 ml volumetric flask, dissolve it with an appropriate amount of ultrapure water, cool it down to room temperature, then dilute it to the mark, shake well, and the solution is obtained.

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1.2.3. Test solution (测试溶液):

Test solution: Take approximately 25 mg of Norfloxacin Lactate, accurately weigh it, place it in a 25 ml volumetric flask, add 2 ml of 0.1 M hydrochloric acid solution to dissolve it, dilute with water to the mark, and shake well; accurately measure 1.5 ml of the above solution, place it in a 50 ml volumetric flask, dilute with water to the mark, shake well, and the solution is obtained.

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Reference substance solution: Take approximately 25 mg of norfloxacin reference substance, accurately weigh it, place it in a 25 ml volumetric flask, add 2 ml of 0.1 M hydrochloric acid solution to dissolve it, dilute with water to the mark, and shake well; accurately measure 1.5 ml of the above solution, place it in a 50 ml volumetric flask, dilute with water to the mark, shake well,

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

2. Chromatogram and data (谱图和数据)

(1) Analysis spectrum of blank solvent (空白溶剂分析图谱):

(2) Analysis spectrum of Test solution (供试品溶液分析图谱):

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(3) Analysis spectrum of reference solution (对照溶液分析图谱):

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
3. Conclusion (结论)

When analyzing the content of norfloxacin lactate under this chromatographic condition using the Welch Ultimate XB-C18 chromatographic column (4.6×150mm, 3μm), the fronting issue of the target compound peak can be effectively improved. This meets the customer's testing requirements.

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