Test report

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Sample Information 样品信息			
Sample Name 样品名称	Styrax 苏合香	Sample properties 样品性状	Thick liquid 浓稠液体
Date of receipt 收样日期	2025/08/12	During testing 测试期间	2025/09/09 ~ 09/10
Test the ingredients and structural formula 测试成分及结构式			
Cinnamic acid 肉桂酸			
Experimental Requirements 实验要求			
The theoretical plate count calculated according to the cinnamic acid peak should not be less than 7000, and the separation degree meets the requirements 理论板数按肉桂酸峰计算应不低于 7000, 分离度满足要求			
Reference Method 参考方法			
China Pharmacopoeia 2025 中国药典 2025			
Reagent information 试剂信息			
Reagent Name 试剂名称	Grade 级别	Brand 品牌	
Glacial acetic acid 冰醋酸	AR	Traditional Chinese medicine 国药	

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Web: www.welchmat.com

1.2.1. Preparation of mobile phase (流动相配制):

Take 500mL methanol and add 500mL 1% acetic acid solution, mix well and ultrasonic degassing is obtained;

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1.2.2. Contrast solution (对照溶液):

Take an appropriate amount of cinnamon acid reference substance, weigh it accurately, and add methanol to make a solution containing 25 μ g per 1ml.

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
1.2.3. Test sample solution (供试品溶液):

Take approximately 0.5g of the sample and add 10ml of freshly prepared potassium hydroxide ethanol solution. Heat and reflux for 1 hour, then rapidly evaporate the ethanol at low temperature. Add 20ml of hot water to the residue for uniform dispersion, let it cool, then add 30ml of water and 20ml of magnesium sulfate solution (1.5 \rightarrow 50). Mix thoroughly, let stand for 10 minutes, filter, and wash the residue with 20ml of water in batches. Combine the filtrate and wash liquid, add hydrochloric acid to form an acidic solution, then extract four times with 40ml of ether each time. Combine the ether extracts, recover the solvent until dry. Dissolve the residue in methanol, transfer to a 100ml volumetric flask, dilute to mark, and shake well. Precisely measure 1ml, place in a 50ml volumetric flask, add methanol to dilute to mark, shake well, filter, and collect the filtrate to obtain the final product.

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

3. Conclusion (结论)

Using Welch Ultimate® Plus-C18 (4.6×250mm,5µm), under this chromatographic condition, the theoretical plate number calculated by cinnamic acid peak is more than 7000, and the separation degree of each peak can meet the test requirements.

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