
Test report

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Sample Information			
样品信息			
Sample Name 样品名称	Lusutrombopag 芦曲泊帕	Sample properties 样品性状	Liquid 液体
Date of receipt 收样日期	2025/09/28	During testing 测试期间	2025/10/11 ~ 2025/10/17
Test the ingredients and structural formula			
测试成分及结构式			
			
Experimental Requirements			
实验要求			
<p>Develop a method to analyze the chiral isomers in the active pharmaceutical ingredient (API) of Lusutrombopag, with requirements that the resolution is greater than 2 and the reproducibility is good.</p> <p>开发方法, 分析芦曲泊帕原料药中手性异构体, 要求分离度大于 2 且重现性良好</p>			
Reference Method			
参考方法			
/			
Reagent information			
试剂信息			
Reagent Name 试剂名称	Grade 级别	Brand 品牌	
Hexane 正己烷	HPLC 色谱级	Aladdin 阿拉丁	
Ethanol	HPLC	Aladdin	

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## 1.2. Solution preparation (溶液配制)

### 1.2.1. Preparation of mobile phase (流动相配制):

Mobile phase: Pipette n-hexane, ethanol, trifluoroacetic acid, and triethylamine, prepare them in the ratio of 950:50:0.5:0.1, and mix them thoroughly to obtain the desired solution.

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### 1.2.2. Blank solution、Diluent (空白溶液、稀释剂):

Ethanol (乙醇)

### 1.2.3. System Solution: Provided by the customer, which can be directly used for analysis.

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
### 1.2.4. Imp reference solution: Provided by the customer, which can be directly used for analysis.

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Tel: 400-810-6969

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
Web: www.welchmat.com


## (2) Analysis spectrum of SST solution (系统溶液分析图谱):

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## (3) Analysis Multi-injection Overlay spectrum of reference solution (多针叠加图谱):

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
Web:www.welchmat.com

### 3. Conclusion (结论)

Using Welch Ultimate Cellu-D (4.6×250mm, 5μm) analyzed the content of lusutrombopag under this chromatographic condition, the resolution between the target compound and its chiral isomers was 2.02, and the results of multiple injections were stable, meeting the customer's testing requirements.

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