Test report

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Sample Information 样品信息			
Sample Name 样品名称	Amoxicillin Sodium and Sulbactam Sodium for Injection 注射用氨苄西林钠舒巴坦钠	Sample properties 样品性状	White solid powder 白色固体粉末
Date of receipt 收样日期	2025/10/20	During testing 测试期间	2025/10/27~10/28
Test the ingredients and structural formula 测试成分及结构式			
/			
Experimental Requirements 实验要求			
Select a C18 chromatographic column that meets the content determination requirements, with a tailing factor between 0.8 and 1.8. 筛选 C18 色谱柱, 满足含量测定要求, 拖尾因子在 0.8-1.8 之间			
Reference Method 参考方法			
Customer Method 客户方法			
Reagent information 试剂信息			
Reagent Name 试剂名称	Grade 级别	Brand 品牌	
Acetonitrile 乙腈	HPLC 色谱纯	Welch 月旭	

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Web: www.welchmat.com

Report No (报告编号): 20250929-1020-03

0.02mol/L Sodium dihydrogen phosphate (pH 4.0): Take 2.4g of sodium dihydrogen phosphate, add it to 1000mL of ultrapure water, dissolve it, filter it, adjust the pH to 4.0 with phosphoric acid, and degas it with ultrasound.

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Acetonitrile: Take an appropriate amount of chromatographically pure acetonitrile and sonicate to obtain it.

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1.2.2. Single standard solution (单标溶液)

Ampicillin solution: Take an appropriate amount of ampicillin and mix it with a flow to form a 1mg/ml solution.

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Sulbactam solution: Take an appropriate amount of Sulbactam and mix it with a flow to form a 1mg/mL solution.

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1.2.3. Reference solution (对照品溶液)

Take the above two solutions and prepare a mixed solution containing 0.6mg/ml ampicillin and 0.3mg/ml sulbactam.

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1.2.4. System suitability solution (系统适用性溶液)

Take 6mg of ampicillin reference substance and 3mg of sulbactam reference substance, dissolve them in 10mL of 0.01mol/L sodium hydroxide, and let them stand at room temperature for 30 minutes. Adjust the pH to 4.0 with phosphoric acid solution;

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Take 5mL of each of the above two solutions and place them in a 25mL volumetric flask containing 5mg of ampicillin and 2.5mg of sulbactam. Dissolve and dilute with mobile phase to the mark.

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Web: www.welchmat.com

1.2.5. Test solution (供试品溶液)

Take an appropriate amount of this product, dissolve it in mobile phase and dilute it to make a solution containing approximately 0.6mg ampicillin and 0.3mg sulbactam per 1mL, filter it, and it is ready.

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