Test report

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Sample Information			
样品信息			
Sample Name 样品名称	lecithin 卵磷脂	Sample properties 样品性状	Liquid 液体
Date of receipt 收样日期	2025/09/28	During testing 测试期间	2025/10/13 ~ 2025/10/15
Test ingredients and structural formula			
测试成分及结构式			
Experimental Requirements			
实验要求			
Develop method to determine egg yolk phosphatidylcholine in phospholipids, requiring that the reference substance does not split peak. 开发方法测定卵磷脂中蛋黄磷脂酰胆碱, 要求对照品不裂峰			
Reference Method			
参考方法			
Method for determination of egg yolk lecithin content in Part IV of the Chinese Pharmacopoeia 中国药典四部蛋黄卵磷脂含量测定方法			
Reagent information			
试剂信息			
Reagent Name 试剂名称	Grade 级别	Brand 品牌	
Water 水	Ultrapure water 超纯水	Welch 月旭	
Methanol 甲醇	HPLC 色谱级	Welch 月旭	

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1.2. Solution preparation (溶液配制)

1.2.1. Preparation of mobile phase (流动相配制):

Mobile phase A: Mix 850ml of methanol and 150ml of water, then add 450 μ l of glacial acetic acid and 50 μ l of triethylamine, mix well, and degas with ultrasound.

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Mobile phase B: Take 100ml of n-hexane, 240ml of isopropanol, and 60ml of mobile phase A 160ml, mix well, and degas with ultrasound.

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

Web: www.welchmat.com

3. Conclusion (结论)

Using the Welch chromatographic column Ultimate SiO₂ (4.6 × 250mm, 5 μ m), analyze the target substance under these chromatographic conditions. The results meet the detection requirements.

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