Test Report

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Sample Information 样品信息			
Sample Name 样品名称	Ampicillin sodium 氨苄西林钠	Sample properties 样品性状	White powder 白色粉末
Date of receipt 收样日期	2025/10/11	During testing 测试期间	2025/10/23~10/24
Test the ingredients and structural formula 测试成分及结构式			
/			
Experimental Requirements 实验要求			
Screen the chromatographic column, with the main peak tailing factor required to be between 0.8 and 1.8. 筛选色谱柱, 要求主峰拖尾因子在 0.8-1.8 之间。			
Reference Method 参考方法			
Customer Method 客户方法			
Reagent information 试剂信息			
Reagent Name 试剂名称	Grade 级别	Brand 品牌	
Acetonitrile 乙腈	HPLC 色谱纯	Welch 月旭	
Potassium dihydrogen	AR	Sinopharm	

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## 1.2. Solution preparation (溶液配制)

### 1.2.1. Preparation of mobile phase (流动相配制)

Mobile phase A: Take 2.72g of potassium dihydrogen phosphate, add it into 100mL of ultrapure water, mix well, and filter. Take 50mL of the above solution and 50mL of chromatographic grade acetonitrile to mix, and then add 900mL of ultrapure water, shake well, and add 0.5mL of 12% acetic acid solution, shake well, and degas by ultrasonic waves.

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Mobile phase B: Take 2.72g of potassium dihydrogen phosphate, add it into 100mL of ultrapure water, mix well, and filter. Take 50mL of the above solution and 400mL of chromatographic grade acetonitrile to mix, and then add 550mL of ultrapure water, shake well, and add 0.5mL of 12% acetic acid solution, shake well, and degas by ultrasonic waves.

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### 1.2.2. System suitability solution (系统适用性溶液):

Take an appropriate of ampicillin and cefradine, add the mobile phase A, and mix well to prepare a solution containing about 0.3mg of ampicillin and 0.02mg of cefradine per 1mL.

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### 1.2.3. Reference Solution (对照品溶液):

Take an appropriate of ampicillin reference substance, add the mobile phase A, and mix well to prepare a solution containing about 1mg of ampicillin per 1mL.

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
### 1.2.4. Test solution (供试品溶液):

Take an appropriate of ampicillin sodium test sample, add the mobile phase A, and mix well to prepare a solution containing about 1mg of ampicillin per 1mL.

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
## 2. Chromatogram and data (谱图和数据)

### (1) System suitability solution (系统适用性溶液)

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
### 3. Conclusion (结论)

Using Welch Ultimate® Polar RP (4.6×250mm, 5μm), the detection requirements of the customer can be met under these chromatographic conditions.

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