Test report

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Sample Information 样品信息			
Sample Name 样品名称	Pioglitazone metformin tablets 吡格列酮二甲双胍片	Sample properties 样品性状	Solid/Liquid 固体/液体
Date of receipt 收样日期	2025/11/17	During testing 测试期间	2025/11/17~11/24
Test the ingredients and structural formula 测试成分及结构式			
/			
Experimental Requirements 实验要求			
System adaptability requirements: The theoretical plate count, calculated based on the pioglitazone peak, should be no less than 10,000, and calculated based on metformin, no less than 6,000. In the reference substance solution, the tailing factors of the pioglitazone peak and the metformin peak were not greater than 2.0 系统适应性要求: 理论板数按吡格列酮峰计算不低于 10000, 按二甲双胍计算不低于 6000; 对照品溶液中, 吡格列酮峰和二甲双胍峰的拖尾因子不大于 2.0			
Reference Method 参考方法			
Customer 客户			
Reagent information 试剂信息			
Reagent Name 试剂名称	Grade 级别	Brand 品牌	Item Number 货号

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1.2. Solution preparation (溶液配制)

1.2.1. Preparation of mobile phase (流动相配制)

Mobile phase: Take 5.75g of ammonium dihydrogen phosphate, add it to 1000mL of ultrapure water, mix well, filter by suction, and degassing by ultrasonic waves; Take 500mL of the above solution and mix it with 500mL of chromatographic grade acetonitrile, shake well, then add 7.2g of sodium dodecyl sulfate, mix well, and ultrasonically degas to obtain the product.

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0.1N hydrochloric acid solution/methanol = 1/1.

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1.2.3. Blank solution/dissolution medium (空白溶液/溶出介质):

Take 19g of citric acid and 2.58g of anhydrous disodium hydrogen phosphate, add them to 1000mL of ultrapure water, mix well, filter by suction, and degas by ultrasonic waves, and the product is obtained.

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1.2.4. Reference substance solution (对照品溶液):

Take 4.6mg of pioglitazone hydrochloride reference substance, place it in a 10mL volumetric flask, add the diluent and dissolve it to the mark, then shake well. Measure 0.4mL and place it in a 10mL volumetric flask. Add 5.6mg of metformin hydrochloride reference substance, dissolve it in the dissolution medium and dilute to the mark, shake well, and it is ready.

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
1.2.5. Test solution (供试品溶液):

Take one tablet of this product, and use 150 ml of pH2.5 buffer solution as the dissolution medium. Ultrasonically treat for 30 minutes until the sample dissolves. Take the solution and continue to filtrate to obtain the product.

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

(3) Diagram of the test solution (供试品溶液图)

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Web:www.welchmat.com

3. Conclusion (结论)

Using Welch Xtimate® UHPLC C18 (2.1×100mm, 1.8μm), the detection requirements of the customer can be met under these chromatographic conditions.

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