

Test report

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Sample Information			
样品信息			
Sample Name 样品名称	Adipic Acid, Trans-2-Butene-1,4-dicarboxylic Acid, 3-Hexenedioic Acid 己二酸, 反-2-丁烯-1,4-二甲酸, 3-己烯二酸	Sample properties 样品性状	Solid 固体
Date of receipt 收样日期	2025/09/29	During testing 测试期间	2025/10/10 ~ 2025/10/15
Test the ingredients and structural formula			
测试成分及结构式			
			
Adipic Acid 己二酸	Trans-2-Butene-1,4-dicarboxylic Acid 反-2-丁烯-1,4-二甲酸	3-Hexenedioic Acid 3-己烯二酸	
Experimental Requirements			
实验要求			
Development of a Method for analyzing adipic acid, trans-2-Butene-1,4-dicarboxylic Acid, and 3-hexenedioic acid, with the requirement of mutual non-Interference and a resolution Greater than 1.5 开发方法分析己二酸, 反-2-丁烯-1,4-二甲酸, 3-己烯二酸, 要求互不干扰分离度大于 1.5			
Reference Method			
参考方法			
/			
Reagent information			
试剂信息			
Reagent Name 试剂名称	Grade 级别	Brand 品牌	

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1.2. Solution preparation (溶液配制)

1.2.1. Preparation of mobile phase (流动相配制):

Mobile phase : Measure 0.5 ml of sulfuric acid into 900 ml of purified water, mix thoroughly.

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0.1M KOH solution

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1.2.3. Test solution (供试品溶液):

Sample solution ,Sample dilution solution and Blank solution are provided by the customer.

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1.2.4. Adipic acid Reference Solution: Place 10 mg of propylene glycol into a 10 mL volumetric flask, dissolve with water, and dilute to the mark.

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1.2.5. Trans-2-Butene-1,4-dicarboxylic acid reference solution: Place 10 mg of propylene glycol into a 10 mL volumetric flask, dissolve with water, and dilute to the mark.

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
1.2.6. 3-Hexenedioic acid reference solution: Place 10 mg of the reference substance into a 10 mL volumetric flask, dissolve with water, and dilute to the mark.

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Tel:400-810-6969

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
Web:www.welchmat.com


(2) Analysis spectrum of Adipic acid solution (己二酸对照溶液分析图谱):

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(3) Analysis spectrum of Trans-2-Butene-1,4-dicarboxylic acid reference solution

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(4) Analysis spectrum of 3-Hexenedioic acid reference solution

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

(5) Analysis spectrum of the Sample solution (供试品溶液分析图谱):

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(6) Analysis spectrum of the Sample dilution solution (供试品稀释溶液分析图谱):

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
3. Conclusion (结论)

Using Welch Xtimate Sugar-H (7.8×300mm, 5µm) analyzed the three target analytes under this chromatographic condition. Trans-2-Butene-1,4-dicarboxylic acid, adipic acid and 3-hexenedioic acid eluted sequentially with the resolution between consecutive peaks being 2.0 and 2.5 respectively, meets the customer's testing requirements.

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